

# Tutorial on how to run “impedimetric\_analysis.py”

## Prerequisites




- Python 3.7.0 or a newer version
- Following packages: *Matplotlib, Numpy, Openpyxl, Os, Pandas, Pylab, Re, Scikit-learn, Scipy, Sys, Tkinter (tk)*.

It is highly recommended that programming beginners should use Anaconda<sup>1</sup>, in which e.g. Spyder can be used to run Python. Anaconda is easy to install and include various packages and tools.

A package can be installed in Anaconda, follow this tutorial<sup>2</sup>. If using Anaconda with Python 3.7.0 or a newer version, following packages are already installed: *Matplotlib, Numpy, Os, Pandas, Re, Scikit-learn Scipy* and *sys*.





To install remaining packages, open an Anaconda prompt/terminal and type following:







- `conda install openpyxl`
- `conda install pylab`
- `conda install tk`

How to run “impedimetric_analysis.py”		
Step no.	What to do	Image
1	<ul style="list-style-type: none"><li>• Open Github and access “Impedimetric_analysis” repository</li><li>• Click “Code” and select “Download ZIP”</li></ul>	
2	<ul style="list-style-type: none"><li>• Unzip the “ZIP” file and open “impedimetric_analysis.py” e.g. via Spyder</li></ul>	
3	<ul style="list-style-type: none"><li>• Run the code</li></ul>	

<sup>1</sup> <https://www.anaconda.com>

<sup>2</sup> <https://docs.anaconda.com/anaconda/user-guide/tasks/install-packages/>

4	<ul style="list-style-type: none"> <li>A pop-up appears</li> <li>Select the analysis method to use</li> </ul> <p>NOTE: For this tutorial, EIS-analysis is selected, but the procedure is the same for the other techniques.</p>	
5	<ul style="list-style-type: none"> <li>Change the settings if relevant</li> <li>Once all settings are set, click “Run!”</li> </ul>	
6	<ul style="list-style-type: none"> <li>A pop-up appears</li> <li>Select all the excel/csv files to be analyzed by this technique</li> </ul> <p>Note, if EIS-analysis is chosen, choose EIS-measurement files only.</p>	
7	<ul style="list-style-type: none"> <li>When the files are selected, click “Open”</li> </ul>	



8	<ul style="list-style-type: none"> <li>Pop-up text “Done” appears</li> <li>Click “OK”</li> </ul>	
9	<ul style="list-style-type: none"> <li>Pop-up text “Choose saving folder” appears</li> <li>Click “OK”</li> </ul>	
10	<ul style="list-style-type: none"> <li>Pop-up appears</li> <li>Select the folder for which the output plots shall be saved</li> <li>Once selected, click “Choose”</li> </ul>	
11	<ul style="list-style-type: none"> <li>Pop-up text “Done” appears</li> <li>Click “OK”</li> </ul>	
12	<ul style="list-style-type: none"> <li>You can now close the initial pop-up by clicking the red cross in the upper left</li> </ul>	
13	<ul style="list-style-type: none"> <li>The plots can be found in your selected folder</li> </ul>	

## Comments on the different plots obtained from EIS-analysis

A total of four different plots are obtained in a successful EIS analysis.

See an example below.

See an example below.

Plot	Comment																									
<div><p>Overview of Resist Charge Transfer for each file</p></div>	A bar chart of the estimated semicircle diameter for each of the selected files.																									
<div><table><thead><tr><th></th><th>Resist Charge Transfer / <math>\Omega</math></th></tr></thead><tbody><tr><td>BSA</td><td>1530.0</td></tr><tr><td>TSH</td><td>1580.0</td></tr><tr><td>TSH-antibody</td><td>686.0</td></tr><tr><td>Clean</td><td>42.0</td></tr></tbody></table></div>		Resist Charge Transfer / $\Omega$	BSA	1530.0	TSH	1580.0	TSH-antibody	686.0	Clean	42.0	A table of the estimated semicircle diameter for each of the selected files.															
	Resist Charge Transfer / $\Omega$																									
BSA	1530.0																									
TSH	1580.0																									
TSH-antibody	686.0																									
Clean	42.0																									
<div><table><thead><tr><th></th><th>BSA Resist Charge Transfer: 1530.0 <math>\Omega</math></th><th>TSH Resist Charge Transfer: 1580.0 <math>\Omega</math></th><th>TSH-antibody Resist Charge Transfer: 686.0 <math>\Omega</math></th><th>Clean Resist Charge Transfer: 42.0 <math>\Omega</math></th></tr></thead><tbody><tr><td>BSA Resist Charge Transfer: 1530.0 <math>\Omega</math></td><td>0.0</td><td>50.0</td><td>-844.0</td><td>-1488.0</td></tr><tr><td>TSH Resist Charge Transfer: 1580.0 <math>\Omega</math></td><td>-50.0</td><td>0.0</td><td>-894.0</td><td>-1538.0</td></tr><tr><td>TSH-antibody Resist Charge Transfer: 686.0 <math>\Omega</math></td><td>844.0</td><td>894.0</td><td>0.0</td><td>-644.0</td></tr><tr><td>Clean Resist Charge Transfer: 42.0 <math>\Omega</math></td><td>1488.0</td><td>1538.0</td><td>644.0</td><td>0.0</td></tr></tbody></table></div>		BSA Resist Charge Transfer: 1530.0 $\Omega$	TSH Resist Charge Transfer: 1580.0 $\Omega$	TSH-antibody Resist Charge Transfer: 686.0 $\Omega$	Clean Resist Charge Transfer: 42.0 $\Omega$	BSA Resist Charge Transfer: 1530.0 $\Omega$	0.0	50.0	-844.0	-1488.0	TSH Resist Charge Transfer: 1580.0 $\Omega$	-50.0	0.0	-894.0	-1538.0	TSH-antibody Resist Charge Transfer: 686.0 $\Omega$	844.0	894.0	0.0	-644.0	Clean Resist Charge Transfer: 42.0 $\Omega$	1488.0	1538.0	644.0	0.0	A table that compares the semicircle diameter for all the files and show the differences.
	BSA Resist Charge Transfer: 1530.0 $\Omega$	TSH Resist Charge Transfer: 1580.0 $\Omega$	TSH-antibody Resist Charge Transfer: 686.0 $\Omega$	Clean Resist Charge Transfer: 42.0 $\Omega$																						
BSA Resist Charge Transfer: 1530.0 $\Omega$	0.0	50.0	-844.0	-1488.0																						
TSH Resist Charge Transfer: 1580.0 $\Omega$	-50.0	0.0	-894.0	-1538.0																						
TSH-antibody Resist Charge Transfer: 686.0 $\Omega$	844.0	894.0	0.0	-644.0																						
Clean Resist Charge Transfer: 42.0 $\Omega$	1488.0	1538.0	644.0	0.0																						
<div><p>Impedance Nyquist</p></div>	A plot of all the selected files.																									